

PowerFlex[®]P3

拡張圧力	バルーン径 (mm)							
	4	5	6	7	8	9	10	12
4atm	-	-	-	-	-	-	-	11.3
6atm	3.8	4.7	5.7	6.6	7.6	8.4	9.4	12.0
8atm	3.9	4.9	5.8	6.8	7.8	8.7	9.7	12.4
10atm	4.0	5.0	6.0	7.0	8.0	9.0	10.0	-
12atm	4.1	5.1	6.1	7.2	8.2	9.2	10.2	-
14atm	4.2	5.2	6.2	7.3	8.4	9.4	10.4	-
15atm	4.2	5.3	6.3	7.4	8.4	-	-	-

Nominal

Rated Burst

■承認番号:21300BZY00256000 ■販売名:拡張用PTAカテーテル
 ■製造販売元: Cardinal Health Japan合同会社

Cordis[®]
 A Cardinal Health company