

造影用カテーテル

■ Tempo[®] 4・Tempo[®] 5 テンポ[®] 4・テンポ[®] 5


















■ 特定保険医療材料請求区分：血管造影用カテーテル／一般用

■ 保険償還価格：¥2,000 / 1本

■ クラス分類 IV

■ JMDN コード：10688104

□ 推奨ガイドワイヤー径：0.038"

(外径/内径) 4F / 0.042"	(外径/内径) 5F / 0.048"	有効長 (cm)	先端形状		side holes
カタログ番号	カタログ番号				
451-413H0	451-513H0	100	BERN		0
451-414H0	451-514H0	100	VER		0
—	451-523H0	100	JB1		0
451-424H0	451-524H0	100	JB2		0
451-430H0	451-530H0	100	SIM1		0
—	451-530H2**	100	SIM1		2
451-431H0	451-531H0**	100	SIM2		0
451-432H0	—	100	SIM3		0
451-435H0	451-535H0	100	H1		0
451-442V0	—	65	C1		0
451-443V0	451-543V0	65	C2		0
451-443H0	451-543H0	100	C2		0
451-446S0	451-546S0	80	RDC (Small)		0
451-447S0	451-547S0	80	RDC		0
451-449V0	451-549V0	65	SHK0.8		0
451-450V0	451-550V0	65	SHK1.0		0
451-456V0	451-556V0	65	J2		0

備考：5本入／箱

承認番号：21000BZY00228000 販売名：コーディス 血管造影用カテーテル

**印の製品は在庫終了次第、販売終了となります。

■ Tempo[®] 4・Tempo[®] 5 テンポ[®] 4・テンポ[®] 5











■特定保険医療材料請求区分：血管造影用カテーテル／一般用

■保険償還価格：¥2,000／1本

■クラス分類 IV

■JMDNコード：10688104

□推奨ガイドワイヤー径：0.038"

(外径/内径) 4F / 0.042"	(外径/内径) 5F / 0.048"	有効長 (cm)	先端形状		side holes
カタログ番号	カタログ番号				
451-457V0	451-557V0	65	J3		0
SRD-5497†	—	80	RH		0
SRD-5821	—	65	VER		0
SRD5627	—	125	VER		0
SRD5498	—	65	STR		0
SRD-5641**	—	65	SHK MOD		0
SRD-6012	SRD-6013	125	JB2		0
SRD6885	—	40	STR		0
SRD6886	—	40	JR4		0
SRD6887	—	40	BER2		0










備考：5本入／箱

承認番号：21000BZY00228000 販売名：コーディス 血管造影用カテーテル

†SRD-5497のみ3本入／箱となります。

**印の製品は在庫終了次第、販売終了となります。

□推奨ガイドワイヤー径：0.035"

(外径/内径) 4F / 0.042"	(外径/内径) 5F / 0.048"	有効長 (cm)	先端形状		side holes
カタログ番号	カタログ番号				
451-401V5	451-501V5	65	STR		5
451-401F5	—	90	STR		5
451-401H5	451-501H5	100	STR		5
451-401L5	—	110	STR		5
451-403V5	451-503V5	65	PIG		5
451-403F5	451-503F5	90	PIG		5
451-403L5	451-503L5	110	PIG		5
451-404V5	451-504V5	65	UNIV		5
451-404F5	—	90	UNIV		5

備考：5本入／箱

承認番号：21000BZY00228000 販売名：コーディス 血管造影用カテーテル