

# THRILL SLALOM™

拡張圧力	バルーン径 (mm)					
	3	4	5	6	7	8
4atm	-	-	-	-	6.5	7.7
6atm	2.8	3.8	4.8	5.7	6.8	8.0
8atm	2.9	3.9	4.9	5.9	7.0	8.3
10atm	3.0	4.0	5.0	6.0	7.3	8.6
12atm	3.1	4.1	5.1	6.1	7.5	-
14atm	3.2	4.2	5.2	6.3	-	-
16atm	3.2	4.4	5.4	-	-	-
18atm	3.3	4.5	-	-	-	-
20atm	3.4	-	-	-	-	-

**Nominal**

**Rated Burst**

■承認番号: 21300BZY00203000 ■販売名: スロラーム  
 ■製造販売元: Cardinal Health Japan 合同会社

**Cordis**  
 A Cardinal Health company